

Sample: \_\_\_\_\_ Sample Concentration: \_\_\_\_\_  
 Sample Buffer: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reservoir Volume: \_\_\_\_\_ Temperature: \_\_\_\_\_  
 Drop Volume: Total \_\_\_\_\_ µl Sample \_\_\_\_\_ µl Reservoir \_\_\_\_\_ µl Additive \_\_\_\_\_ µl

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

### Crystal Screen HT™ - HR2-130 Scoring Sheet

Date: \_\_\_\_\_ Date: \_\_\_\_\_ Date: \_\_\_\_\_

1. (A1)	0.02 M Calcium chloride dihydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
2. (A2)	0.4 M Potassium sodium tartrate tetrahydrate			
3. (A3)	0.4 M Ammonium phosphate monobasic			
4. (A4)	0.1 M TRIS hydrochloride pH 8.5, 2.0 M Ammonium sulfate			
5. (A5)	0.2 M Sodium citrate tribasic dihydrate, 0.1 M HEPES sodium pH 7.5, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
6. (A6)	0.2 M Magnesium chloride hexahydrate, 0.1 M TRIS hydrochloride pH 8.5, 30% w/v Polyethylene glycol 4,000			
7. (A7)	0.1 M Sodium cacodylate trihydrate pH 6.5, 1.4 M Sodium acetate trihydrate			
8. (A8)	0.2 M Sodium citrate tribasic dihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 30% v/v 2-Propanol			
9. (A9)	0.2 M Ammonium acetate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 30% w/v Polyethylene glycol 4,000			
10. (A10)	0.2 M Ammonium acetate, 0.1 M Sodium acetate trihydrate pH 4.6, 30% w/v Polyethylene glycol 4,000			
11. (A11)	0.1 M Sodium citrate tribasic dihydrate pH 5.6, 1.0 M Ammonium phosphate monobasic			
12. (A12)	0.2 M Magnesium chloride hexahydrate, 0.1 M HEPES sodium pH 7.5, 30% v/v 2-Propanol			
13. (B1)	0.2 M Sodium citrate tribasic dihydrate, 0.1 M TRIS hydrochloride pH 8.5, 30% v/v Polyethylene glycol 400			
14. (B2)	0.2 M Calcium chloride dihydrate, 0.1 M HEPES sodium pH 7.5, 28% v/v Polyethylene glycol 400			
15. (B3)	0.2 M Ammonium sulfate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 30% w/v Polyethylene glycol 8,000			
16. (B4)	0.1 M HEPES sodium pH 7.5, 1.5 M Lithium sulfate monohydrate			
17. (B5)	0.2 M Lithium sulfate monohydrate, 0.1 M TRIS hydrochloride pH 8.5, 30% w/v Polyethylene glycol 4,000			
18. (B6)	0.2 M Magnesium acetate tetrahydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 20% w/v Polyethylene glycol 8,000			
19. (B7)	0.2 M Ammonium acetate, 0.1 M TRIS hydrochloride pH 8.5, 30% v/v 2-Propanol			
20. (B8)	0.2 M Ammonium sulfate, 0.1 M Sodium acetate trihydrate pH 4.6, 25% w/v Polyethylene glycol 4,000			
21. (B9)	0.2 M Magnesium acetate tetrahydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
22. (B10)	0.2 M Sodium acetate trihydrate, 0.1 M TRIS hydrochloride pH 8.5, 30% w/v Polyethylene glycol 4,000			
23. (B11)	0.2 M Magnesium chloride hexahydrate, 0.1 M HEPES sodium pH 7.5, 30% v/v Polyethylene glycol 400			
24. (B12)	0.2 M Calcium chloride dihydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 20% v/v 2-Propanol			
25. (C1)	0.1 M Imidazole pH 6.5, 1.0 M Sodium acetate trihydrate			
26. (C2)	0.2 M Ammonium acetate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
27. (C3)	0.2 M Sodium citrate tribasic dihydrate, 0.1 M HEPES sodium pH 7.5, 20% v/v 2-Propanol			
28. (C4)	0.2 M Sodium acetate trihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 30% w/v Polyethylene glycol 8,000			
29. (C5)	0.1 M HEPES sodium pH 7.5, 0.8 M Potassium sodium tartrate tetrahydrate			
30. (C6)	0.2 M Ammonium sulfate, 30% w/v Polyethylene glycol 8,000			
31. (C7)	0.2 M Ammonium sulfate, 30% w/v Polyethylene glycol 4,000			
32. (C8)	2.0 M Ammonium sulfate			
33. (C9)	4.0 M Sodium formate			
34. (C10)	0.1 M Sodium acetate trihydrate pH 4.6, 2.0 M Sodium formate			
35. (C11)	0.1 M HEPES sodium pH 7.5, 0.8 M Sodium phosphate monobasic monohydrate, 0.8 M Potassium phosphate monobasic			
36. (C12)	0.1 M TRIS hydrochloride pH 8.5, 8% w/v Polyethylene glycol 8,000			
37. (D1)	0.1 M Sodium acetate trihydrate pH 4.6, 8% w/v Polyethylene glycol 4,000			
38. (D2)	0.1 M HEPES sodium pH 7.5, 1.4 M Sodium citrate tribasic dihydrate			
39. (D3)	0.1 M HEPES sodium pH 7.5, 2% v/v Polyethylene glycol 400, 2.0 M Ammonium sulfate			
40. (D4)	0.1 M Sodium citrate tribasic dihydrate pH 5.6, 20% v/v 2-Propanol, 20% w/v Polyethylene glycol 4,000			
41. (D5)	0.1 M HEPES sodium pH 7.5, 10% v/v 2-Propanol, 20% w/v Polyethylene glycol 4,000			
42. (D6)	0.05 M Potassium phosphate monobasic, 20% w/v Polyethylene glycol 8,000			
43. (D7)	30% w/v Polyethylene glycol 1,500			
44. (D8)	0.2 M Magnesium formate dihydrate			
45. (D9)	0.2 M Zinc acetate dihydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 18% w/v Polyethylene glycol 8,000			
46. (D10)	0.2 M Calcium acetate hydrate, 0.1 M Sodium cacodylate trihydrate pH 6.5, 18% w/v Polyethylene glycol 8,000			
47. (D11)	0.1 M Sodium acetate trihydrate pH 4.6, 2.0 M Ammonium sulfate			
48. (D12)	0.1 M TRIS hydrochloride pH 8.5, 2.0 M Ammonium phosphate monobasic			

**HAMPTON**  
**RESEARCH**  
*Solutions for Crystal Growth*

34 Journey  
 Also Viejo, CA 92656-5317 U.S.A.  
 Tel: (949) 425-1321 • Fax: (949) 425-1611  
 e-mail: tech@hampont.com  
 Website: www.hampontresearch.com

Sample: \_\_\_\_\_ Sample Concentration: \_\_\_\_\_  
 Sample Buffer: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reservoir Volume: \_\_\_\_\_ Temperature: \_\_\_\_\_  
 Drop Volume: Total \_\_\_\_\_ µl Sample \_\_\_\_\_ µl Reservoir \_\_\_\_\_ µl Additive \_\_\_\_\_ µl

- 1 Clear Drop
- 2 Phase Separation
- 3 Regular Granular Precipitate
- 4 Birefringent Precipitate or Microcrystals
- 5 Posettes or Spherulites
- 6 Needles (1D Growth)
- 7 Plates (2D Growth)
- 8 Single Crystals (3D Growth < 0.2 mm)
- 9 Single Crystals (3D Growth > 0.2 mm)

### Crystal Screen HT™ - HR2-130 Scoring Sheet

Date:      Date:      Date:

49. (E1)	2.0 M Sodium chloride, 10% w/v Polyethylene glycol 6,000			
50. (E2)	0.5 M Sodium chloride, 0.01 M Magnesium chloride hexahydrate, 0.01 M Hexadecyltrimethylammonium bromide			
51. (E3)	25% v/v Ethylene glycol			
52. (E4)	35% v/v 1,4-Dioxane			
53. (E5)	2.0 M Ammonium sulfate, 5% v/v 2-Propanol			
54. (E6)	1.0 M Imidazole pH 7.0			
55. (E7)	10% w/v Polyethylene glycol 1,000, 10% w/v Polyethylene glycol 8,000			
56. (E8)	1.5 M Sodium chloride, 10% v/v Ethanol			
57. (E9)	0.1 M Sodium acetate trihydrate pH 4.6, 2.0 M Sodium chloride			
58. (E10)	0.2 M Sodium chloride, 0.1 M Sodium acetate trihydrate pH 4.6, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
59. (E11)	0.01 M Cobalt(II) chloride hexahydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 1.0 M 1,6-Hexanediol			
60. (E12)	0.1 M Cadmium chloride hydrate, 0.1 M Sodium acetate trihydrate pH 4.6, 30% v/v Polyethylene glycol 400			
61. (F1)	0.2 M Ammonium sulfate, 0.1 M Sodium acetate trihydrate pH 4.6, 30% w/v Polyethylene glycol monomethyl ether 2,000			
62. (F2)	0.2 M Potassium sodium tartrate tetrahydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 2.0 M Ammonium sulfate			
63. (F3)	0.5 M Ammonium sulfate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 1.0 M Lithium sulfate monohydrate			
64. (F4)	0.5 M Sodium chloride, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 2% v/v Ethylene imine polymer			
65. (F5)	0.1 M Sodium citrate tribasic dihydrate pH 5.6, 35% v/v tert-Butanol			
66. (F6)	0.01 M Iron(III) chloride hexahydrate, 0.1 M Sodium citrate tribasic dihydrate pH 5.6, 10% v/v Jeffamine® M-600®			
67. (F7)	0.1 M Sodium citrate tribasic dihydrate pH 5.6, 2.5 M 1,6-Hexanediol			
68. (F8)	0.1 M MES monohydrate pH 6.5, 1.6 M Magnesium sulfate heptahydrate			
69. (F9)	0.1 M Sodium phosphate monobasic monohydrate, 0.1 M Potassium phosphate monobasic, 0.1 M MES monohydrate pH 6.5, 2.0 M Sodium chloride			
70. (F10)	0.1 M MES monohydrate pH 6.5, 12% w/v Polyethylene glycol 20,000			
71. (F11)	1.6 M Ammonium sulfate, 0.1 M MES monohydrate pH 6.5, 10% v/v 1,4-Dioxane			
72. (F12)	0.05 M Cesium chloride, 0.1 M MES monohydrate pH 6.5, 30% v/v Jeffamine® M-600®			
73. (G1)	0.01 M Cobalt(II) chloride hexahydrate, 0.1 M MES monohydrate pH 6.5, 1.8 M Ammonium sulfate			
74. (G2)	0.2 M Ammonium sulfate, 0.1 M MES monohydrate pH 6.5, 30% w/v Polyethylene glycol monomethyl ether 5,000			
75. (G3)	0.01 M Zinc sulfate heptahydrate, 0.1 M MES monohydrate pH 6.5, 25% v/v Polyethylene glycol monomethyl ether 550			
76. (G4)	1.6 M Sodium citrate tribasic dihydrate pH 6.5			
77. (G5)	0.5 M Ammonium sulfate, 0.1 M HEPES pH 7.5, 30% v/v (+/-)-2-Methyl-2,4-pentanediol			
78. (G6)	0.1 M HEPES pH 7.5, 10% w/v Polyethylene glycol 6,000, 5% v/v (+/-)-2-Methyl-2,4-pentanediol			
79. (G7)	0.1 M HEPES pH 7.5, 20% v/v Jeffamine® M-600®			
80. (G8)	0.1 M Sodium chloride, 0.1 M HEPES pH 7.5, 1.6 M Ammonium sulfate			
81. (G9)	0.1 M HEPES pH 7.5, 2.0 M Ammonium formate			
82. (G10)	0.05 M Cadmium sulfate hydrate, 0.1 M HEPES pH 7.5, 1.0 M Sodium acetate trihydrate			
83. (G11)	0.1 M HEPES pH 7.5, 70% v/v (+/-)-2-Methyl-2,4-pentanediol			
84. (G12)	0.1 M HEPES pH 7.5, 4.3 M Sodium chloride			
85. (H1)	0.1 M HEPES pH 7.5, 10% w/v Polyethylene glycol 8,000, 8% v/v Ethylene glycol			
86. (H2)	0.1 M HEPES pH 7.5, 20% w/v Polyethylene glycol 10,000			
87. (H3)	0.2 M Magnesium chloride hexahydrate, 0.1 M Tris pH 8.5, 3.4 M 1,6-Hexanediol			
88. (H4)	0.1 M Tris pH 8.5, 25% v/v tert-Butanol			
89. (H5)	0.01 M Nickel(II) chloride hexahydrate, 0.1 M Tris pH 8.5, 1.0 M Lithium sulfate monohydrate			
90. (H6)	1.5 M Ammonium sulfate, 0.1 M Tris pH 8.5, 12% v/v Glycerol			
91. (H7)	0.2 M Ammonium phosphate monobasic, 0.1 M Tris pH 8.5, 50% v/v (+/-)-2-Methyl-2,4-pentanediol			
92. (H8)	0.1 M Tris pH 8.5, 20% v/v Ethanol			
93. (H9)	0.01 M Nickel(II) chloride hexahydrate, 0.1 M Tris pH 8.5, 20% w/v Polyethylene glycol monomethyl ether 2,000			
94. (H10)	0.1 M Sodium chloride, 0.1 M BICINE pH 9.0, 20% v/v Polyethylene glycol monomethyl ether 550			
95. (H11)	0.1 M BICINE pH 9.0, 2.0 M Magnesium chloride hexahydrate			
96. (H12)	0.1 M BICINE pH 9.0, 2% v/v 1,4-Dioxane, 10% w/v Polyethylene glycol 20,000			



Solutions for Crystal Growth

34 Journey  
 Aliso Viejo, CA 92656-3317 U.S.A.  
 Tel: (949) 425-1321 • Fax: (949) 425-1611  
 e-mail: tech@hrc.com  
 Website: www.hamptonresearch.com

© 1991-2015 Hampton Research Corp. all rights reserved  
 Printed in the United States of America. This guide or parts thereof may not  
 be reproduced in any form without the written permission of the publishers.